



### Helbing Laboratory eDNA Technical Bulletin

All eDNA tools are validated through a rigorous multi-step evaluation protocol that includes tests of DNA target specificity and amplification sensitivity<sup>1-3</sup>.

#### General eDNA Assay Information

Target Species: Cutthroat Trout (*Oncorhynchus clarkii*)  
Species Code: te-ONCL

eDNA qPCR Tool: eONCL4  
eDNA qPCR Format: TaqMan

Gene Target: MT-ND1  
Published in:

#### eDNA Assay Sensitivity Test Summary using gBlocks™ Synthetic DNA

LOD 0.2 95% CI 0.1-0.4 Copies/Rxn LOQ 0.8 95% CI 0.6-1.4 Copies/Rxn LOB 0 hits/8

LOQ<sub>continuous</sub> 4 Copies/Rxn

Binomial-Poisson model for 8 technical replicates determined using eLowQuant R code<sup>4</sup>. When the LOQ < LOD, use the LOD for the LOQ. Enzyme: QIAcuity

#### eDNA Assay Specificity Test Information

Each qPCR reaction in the specificity assay contained 10 picograms of voucher target gDNA (n=25 technical replicates)

Species	Common Name ( <i>Species</i> )	Detection	# Voucher		Sample Sources/Locations
			Specimens		
ma-HOSA	Human ( <i>Homo Sapiens</i> )	No	1		Netherlands
te-ONCLcl	Coastal Cutthroat Trout ( <i>Oncorhynchus clarkii clarkii</i> )	Yes	5		British Columbia
te-ONCLle	Westslope Cutthroat Trout ( <i>Oncorhynchus clarkii lewisi</i> )	Yes	9		Alberta
te-ONGO	Pink Salmon ( <i>Oncorhynchus gorbuscha</i> )	No	1		British Columbia
te-ONKE	Chum Salmon ( <i>Oncorhynchus keta</i> )	No	1		British Columbia
te-ONKI	Coho Salmon ( <i>Oncorhynchus kisutch</i> )	No	1		British Columbia
te-ONMY	Rainbow Trout ( <i>Oncorhynchus mykiss</i> )	No	6		Alberta and British Columbia
te-ONNE	Sockeye Salmon ( <i>Oncorhynchus nerka</i> )	No	1		British Columbia
te-ONTS	Chinook Salmon ( <i>Oncorhynchus tshawytscha</i> )	No	1		British Columbia
te-SACO	Bull Trout ( <i>Salvelinus confluentus</i> )	No	4		Alberta
te-SAFO	Brook Trout ( <i>Salvelinus fontinalis</i> )	No	4		Alberta
te-SAMA	Dolly Varden ( <i>Salvelinus malma</i> )	No	1		Alberta
te-SASA	Atlantic Salmon ( <i>Salmo salar</i> )	No	1		Nova Scotia

#### References



- Hobbs, J, Adams, IT, Round, JM, Goldberg, CS, Allison, MJ, Bergman, LC, Mirabzadeh, A, Allen, H, Helbing, CC (2020) Revising the range of Rocky Mountain tailed frog, *Ascaphus montanus*, in British Columbia, Canada, using environmental DNA methods. Environmental DNA. 2020; 2: 350-361. <https://doi.org/10.1002/edn3.82>
- Hobbs, J, Round, JM, Allison, MJ, Helbing, CC (2019) Expansion of the known distribution of the coastal tailed frog, *Ascaphus truei*, in British Columbia, Canada, using robust eDNA detection methods. PLOS ONE 14(3): e0213849. <https://doi.org/10.1371/journal.pone.0213849>
- Langlois, VS, Allison, MJ, Bergman, LC, To, TA, and Helbing, CC (2020) The need for robust qPCR-based eDNA detection assays in environmental monitoring and risk assessments. Environmental DNA, 3: 519-527. doi: 10.1002/edn3.164
- Lesperance, M, Allison, MJ, Bergman, LC, Hocking, MD, and Helbing, CC (2021) A statistical model for calibration and computation of detection and quantification limits for low copy number environmental DNA samples. Environmental DNA, 00: 1-12. doi: 10.1002/edn3.220



eDNA Assay Sensitivity Test Details using gBlocks™ synthetic DNA



To calculate tables for different numbers of replicates, raw csv data files can be accessed here:  
<https://onlineacademiccommunity.uvic.ca/helbinglab/edna/>

From 8 Technical Replicates




# Detects	# Copies	SE
0	0.00	0.00
1	0.08	0.08
2	0.16	0.12
3	0.26	0.16
4	0.39	0.22
5	0.55	0.28
6	0.78	0.38
7	1.17	0.58

Determined using eLowQuant R code<sup>4</sup>.



Binomial-Poisson model: no intercept  
Determined using eLowQuant R code<sup>4</sup>.  
Based on a 2 µL DNA input in a total 15 µL reaction

Applied to reactions with 100% positive hits



Efficiency 97%

Field Sample Validation

Known  
Sample Type Presence # Samples Detected Location

Abbreviations			
95% CI	95% Confidence interval	LOQ	Limit of quantification
eDNA	Environmental DNA	MT-ND1	Mitochondrial NADH subunit 1 gene
gDNA	Total genomic DNA extracted from voucher specimen	NTC	qPCR no template control
LOB	Limit of blank	qPCR	Quantitative real-time polymerase chain reaction
LOD	Limit of detection	SE	Standard error